

INTEGRATION TEST DOCUMENT		Application/System Name	NGITIFY: A DENTAL INFORMATION MANAGEMENT SYSTEM WITH AI- DRIVEN PATIENT ENGAGEMENT FOR DENTIME DENTAL CLINIC		
Test Cycle No.	1	User/Access Type	Dentist		
Date Tested	08/25/2026	Type of System	Web Application		
Pre-conditions	Open a supported browser; use prepared accounts and test data for the stated role and branch scope				
Verification Steps	Enter the correct URL, execute every listed action, and record Actual Results and Remarks				
Modules	Action Description		Expected Results	Actual Results	Remarks
Radiograph Upload	Choose an unsupported, oversized, or missing image and click Upload or Save.		A file validation error displays and no radiograph is stored.		
Radiograph Upload	Choose a valid image, enter metadata, and click Cancel.		The upload closes and no new radiograph appears.		
Radiograph Upload	Upload a valid radiograph for the selected patient.		One image is stored under the correct patient and appears in the gallery with its metadata.		
Radiograph Gallery	Open an image and use Original, Enhanced, zoom, Fit, Grab, or Compare controls.		The viewer changes only the presentation or comparison state and preserves the stored original.		
Radiograph Gallery	Click Back or Close in the viewer.		The viewer closes and returns to the selected patient's radiograph gallery.		
Delete Radiograph Modal	Open delete confirmation and click Cancel.		The modal closes and the radiograph remains available.		

Delete Radiograph Modal	Confirm deletion of an eligible radiograph.	The selected image is removed once or a dependency blocker displays, and the action is audited.		
AI Radiograph Review	Request AI enhancement or analysis and wait for completion.	A separate review result or recoverable error displays without replacing the original image or publishing unapproved findings.		
Radiograph Annotation	Add, edit, reset, or remove a visual annotation.	The active annotation layer updates as selected while the original radiograph remains unchanged.		
Dentist Finding Modal	Add or edit a finding and click Cancel.	The modal closes and no finding changes.		
Dentist Finding Modal	Save a valid dentist finding or archive and restore it.	The finding state updates once with dentist attribution and remains linked to the correct image and patient.		
Delete Finding Modal	Click Cancel, then repeat and confirm Delete.	Cancel preserves the finding; Delete removes the selected eligible finding once and refreshes the review.		

Radiograph Review Approval	Approve the completed review summary.	Only dentist-approved findings and explanation become available to the patient; unapproved AI content remains internal.		
Material Usage	Search, filter, and expand material usage logs.	Matching dentist-authorized usage entries display with appointment, item, quantity, and date details.		
Material Usage Modal	Add or remove item rows and click Cancel.	The modal closes and no usage record or stock deduction occurs.		
Material Usage Modal	Submit invalid quantities, insufficient stock, or an ineligible appointment.	Validation displays and inventory remains unchanged.		
Material Usage Modal	Submit valid materials for a completed appointment.	One usage log is created and the correct inventory batches and totals are reduced.		
Dentist Scope	Attempt to open an unassigned patient's EMR through a copied URL.	Access is denied or redirected and no unassigned patient information is exposed.		